



reddot design award
winner 2011

MagniLink Vision



MagniLink Vision

MagniLink Vision is a modular system which allows the user to choose key features. Start by choosing camera, screen size and then what features that meets your preferences. MagniLink Vision makes every day easier, and it also provides endless possibilities. MagniLink Vision is one step ahead...

in design.

With soft edges and a combination of light and dark grey colours, MagniLink Vision can, at first sight, blend in or stand out. A closer inspection reveals that the cords have all been concealed. This creates the harmonious shape. The use of cast aluminium enhances both the quality and the visual impression.

in quality.

One key factor of MagniLink Vision is the high picture quality. The unit offers better resolution, brighter colours and image sharpness that is one step ahead. Our cameras with High Resolution and High Definition are beyond any standard solutions, but they are the standard for MagniLink Vision and provide full digital image transfer.

for users.

Sit down in front of the unit, and start working with your MagniLink Vision. Adjust the height, select the colour and magnification level and enjoy the smooth reading table. For an ideal reading position, the reading table has a rounded front edge that is comfortable to rest your hands on.

To step ahead even further regarding ergonomics, MagniLink Vision has a low weight and offers the possibility to listen to text (MLV-13). All control buttons are integrated in the touchscreen and are easily accessible.



reddot design award
winner 2011



Choose your model

- 1 - Choose your camera
 - Sony High Resolution (HR)
 - Sony High Definition (HD)

- 2 - Choose your screen size
 - 18,5" HD-screen
 - 23" FHD-screen

- 3 - Choose your control panel
 - 3-push Panel: Controls magnification, natural and artificial colours.
 - 2-knob Panel: Has the same features as above but offers a even more tactile experience with knobs.
 - 5-push Panel: Controls more features, such as guiding line, change of contrast, a focus lock and the option to work with a computer.

Further possibilities

- MagniLink Vision-system can be connected to a PC or MAC. As accessory LVI offer a DVI Switch or software that provides different split screen modes and possibility to take a picture or record a video.
- EUC software (End User Configuration) can modify the settings according to the user, such as magnification, guiding line, curtains, artificial colours and lighting.
- MagniLink Vision also offers Text-To-Speech functionality. Simply touch the screen to listen to the text. It will increase your reading speed, accuracy, and endurance.



TECHNICAL SPECIFICATIONS

Model	MagniLink Vision						
	Art nr	MLV-01	MLV-02	MLV-03	MLV-04	MLV-05	MLV-06
Magnification	2,5-50x	2,5-50x	2,5-50x	3,0-50x	3,0-50x	3,0-50x	3,0-50x
Camera	HR	HR	HR	HR	HR	HR	HR
Screen size	18,5"	18,5"	18,5"	23"	23"	23"	23"
Control panel	3-push	2-turn	5-push	3-push	2-turn	5-push	5-push
Guiding line	No	No	Yes	No	No	Yes	Yes
Focus lock	No	No	Yes	No	No	Yes	Yes
Contrast	Auto	Auto	Manual	Auto	Auto	Manual	Manual
Weight	11,8 kg	11,8 kg	11,8 kg	13,5 kg	13,5 kg	13,5 kg	13,5 kg
Text-to-speech	No	No	No	No	No	No	No

Model	MagniLink Vision							
	Art nr	MLV-07	MLV-08	MLV-09	MLV-10	MLV-11	MLV-12	MLV-13
Magnification	1,7-60x	1,7-60x	1,7-60x	2,0-60x	2,0-60x	2,0-60x	2,0-60x	2,0-60x
Camera	HD	HD	HD	HD	HD	HD	HD	HD
Screen size	18,5"	18,5"	18,5"	23"	23"	23"	23"	23"
Control panel	3-push	2-turn	5-push	3-push	2-turn	5-push	5-push	5-push
Guiding line	No	No	Yes	No	No	Yes	Yes	Yes
Focus lock	No	No	Yes	No	No	Yes	Yes	Yes
Contrast	Auto	Auto	Manual	Auto	Auto	Manual	Manual	Manual
Weight	11,8 kg	11,8 kg	11,8 kg	13,5 kg	13,5 kg	13,5 kg	13,5 kg	11,8 kg
Text-to-speech	No	No	No	No	No	No	No	Yes



Eschenbach Optik of America, Inc.
22 Shelter Rock Lane
Danbury, CT 06810
Tel: (800) 487-5389
Fax: (888) 799-7200
www.eschenbach.com
Email: info@eschenbach.com



 Better Vision. Better Life.